

Hei-Standard Hotplate Trade-In Promotion

The MR Hei-Standard Hotplate package provides you with temperature control, strong magnetic field coverage, high temperature capacity, and variable speeds to match your application needs. Trade-In ANY old hotplate for a NEW Heidolph MR Hei-Standard Hotplate system TODAY and start utilizing the benefits of a Heidolph product designed for your Research.



- Kera-Disk coating of 145mm diameter hotplate provides you both quick heating times as observed with standard aluminum plate and chemical resistance to prevent corrosion.
- 600 Watts of heating power gives you uniform heat transfer across your vessel with a rated maximum temperature of 300°C
- An independent safety circuit switches MR Hei-Standard off if unit heats 25°C above set point to ensure your safety is controlled.
- Separate analog knobs control heating or stirring with a separate on/off switch for the heating function to prevent accidental heating incidents
- EKT Hei-Con digital contact thermometer gives you the ability to control your heating via immersion into your vessel or reaction mixture

Trade-In and Upgrade Your Hotplate Today!!!

	Description	List	Customer Cost	Cost Savings
#036111210	MR Hei-Standard Package	\$995.00	\$699.00	\$296.00

* Introductory Offer Expires 12/31/10; Includes MR Hei-Standard Hotplate, EKT Hei-Con, Support Rod, Temperature Sensor Holder.

Contact a local Heidolph Brinkmann Sales rep TODAY to learn more about the Radleys Heat-On Block system that fits directly onto the MR Hei-Standard Hotplate Package

Heidolph Brinkmann, LLC.
1241 Jarvis Avenue Elk Grove Village, IL 60007
Tel 1-866-650-9604 Fax 1-224-265-9611
www.heidolphbrinkmann.com

*"From Life Science to Rocket Science our products
and expertise help support the top scientific minds in America"*