

Reactor-Ready Systems Promotion:

The **Reactor-Ready™ Reaction System** from **Radleys** is a new, innovative and affordable multi-volume jacketed glass reaction station for labs involved with synthesis and process research and development. Simple to use, quick to set-up and incorporating a range of innovative features; the Reactor-Ready is designed to boost a chemist's productivity.



Reactor-Ready™ Reaction System:

Features:

- Wide Range of Jacket Vessel Volumes, 250mL to 5L
- Quick and Easy Exchanging of Vessels without Tools
- Workable Temp Ranges of -60°C to 190°C (dependent on temperature control system)
- Easy to Remove Hose Coupling and Heat Transfer Fluid Manifold Reduces Messy Vessel Exchanges

Complete Systems for 4 Common Vessels Sizes to Make Your Purchasing Simple!!!

Reactor-Ready System Part #	Description	Price
015960100	500mL Jacketed System	\$8,925.00
015960200	1000mL Jacketed System	\$8,995.00
015960300	2000mL Jacketed System	\$9,495.00
015960400	5000mL Jacketed System	\$9,995.00

*All Reactor-Ready Systems include: Reactor-Ready Core, Heidolph RZR 2021 OH Stirrer Jacketed Vessel, Vessel Lid, Stir Shaft, PTFE Temp Probe, and Fittings

**Temperature Control System, Fittings, and Connection Tubing to RR Manifold Sold Separately

*** Pricing Valid till 9/30/2010

Heidolph Brinkmann, LLC.

1241 Jarvis Avenue, Elk Grove Village, IL 60007

Tel 1-866-650-9604 Fax 1-224-265-9611

www.heidolphbrinkmann.com

service@heidolph.com